
 <b>Model Kramp GmbH</b> Otto-Hahn-Straße 41 63456 Hanau Tel. +49 6181-6750-0	Job No.: me4020375	Created at: 07.08.2024	Operator: jb
	Size: 448 x 350 mm (w x h)	1. AC: 11.09.2024 lg	5. AC:
	MN 170002304 MC 0000	2. AC:	6. AC:
	LX 4979 MZ/AZ 1143	3. AC:	7. AC:
	Typestyle: Myriad Pro	Type size: 6 pt	8. AC:
Belotero Revive IfU Russia Russland			

	Released for printing	Procedure no. (or identifier)	Identical to approved texts	Compliant with national blue box requirements	Serialization (readable data)	Non-regulatory information is correct and complete	Date Signature
<b>*By signing off this artwork the reviewer assumes full responsibility for the release of the print. Any malperformance may lead to liability claims.</b>							
Reg. Affairs Reviewer 1 (mandatory)*	<input type="checkbox"/>	...	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> (1)	<input type="checkbox"/> (2)	
Reg. Affairs Reviewer 2 (mandatory)*	<input type="checkbox"/>	...	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> (1)	<input type="checkbox"/> (2)	
Information officer (mandatory)	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		
Distributor (Non-Affiliate countries)	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		
Marketing	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		
Legal for released for printing	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		

(1) e.g., PC-NTIN, lot, expiry date format mm/yyyy; (2) Only for affiliate countries;

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Authorized representative of the manufacturer in the territory of the Russian Federation:  
LLC "Merz Pharma"  
125124, Moscow,  
street 3rd Yamskogo Polya 2,  
bld. 26, 2nd floor, premises VII, room 29;  
Tel.: +7 (495) 653-8555  
fax: +7 (495) 653-8554  
e-mail: info@merz.ru.

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#### Terms of warranty







Anteis S.A., 18 Chemin des Aulx, CH-1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero® Revive guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage and use in accordance with the instructions for use for two years from the date of manufacture.

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	Released for printing	Procedure no. (or identifier)	Identical to approved texts	Compliant with national blue box requirements	Serialization (readable data)	Non-regulatory information is correct and complete	Date Signature
<b>*By signing off this artwork the reviewer assumes full responsibility for the release of the print. Any malperformance may lead to liability claims.</b>							
Reg. Affairs Reviewer 1 (mandatory)*	<input type="checkbox"/>	...	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> (1)	<input type="checkbox"/> (2)	
Reg. Affairs Reviewer 2 (mandatory)*	<input type="checkbox"/>	...	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> (1)	<input type="checkbox"/> (2)	
Information officer (mandatory)	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		
Distributor (Non-Affiliate countries)	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		
Marketing	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		
Legal für released for printing	<input type="checkbox"/> received		<input type="checkbox"/> with comments		<input type="checkbox"/> without comments		

(1) e.g. PC-NTIN, lot, expiry date format mm/yyyy; (2) Only for affiliate countries;

#### GB INSTRUCTIONS FOR USE FOR MEDICAL DEVICE

#### „Viscoelastic implant based on cross-linked sodium hyaluronate Belotero® Revive“

**Product name**  
Viscoelastic implant based on cross-linked sodium hyaluronate Belotero® Revive, consisting of:  
- implantable gel, in a 1.0 ml syringe - 1 pc.;  
- sterile needle 30G ½" (0.30 x 13 mm) (MA No. FSZ 2011/11116 manufactured by TSK Laboratory, Japan) - 2 pcs.;  
- traceability label - 2 pcs.;  
- instructions for use – 1 pc.

#### Intended Use

Belotero® Revive is intended to restore skin quality by revitalizing facial skin through rehydration, improvement of skin elasticity and firmness, and smoothening of superficial fine lines.

#### Intended User

Belotero® Revive is designed to be used in medical and preventive treatment institutions.

#### Target treated population

Adult patients (older than 18 years) regardless of gender who do not meet contraindications.

#### Description

Belotero® Revive is a sterile, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer, containing glycerol. Belotero® Revive is an injectable biodegradable implant. Belotero® Revive degrades gradually in the body. The full degradation of Belotero® Revive following injection (implantation) in the human body is around 6-8 months.

#### Presentation

Belotero® Revive is presented as a single-use pre-filled glass syringe of 1.0 ml containing 1.0 ml of implantable gel, sterilized by moist heat in an autoclave. Belotero® Revive is packaged in a cardboard consumer pack. Each consumer pack includes one instructions for use leaflet and one sealed semi-rigid blister, containing one syringe with sterile implantable gel, two traceability labels and two sterile CE-marked needles for single use.

<b>Composition</b>	
Sodium hyaluronate:	20.0 mg/ml
Glycerol Phosphate Buffer pH 7:	q.s.p. 1.0 ml
<b>Glycerol Phosphate Buffer composition:</b>	
Glycerol (C <sub>3</sub> H <sub>8</sub> O <sub>3</sub> ):	19.0 mg/ml
Disodium hydrogen phosphate dodecahydrate (Na <sub>2</sub> HPO <sub>4</sub> × 12H <sub>2</sub> O):	5.5 mg/ml
Sodium dihydrogen phosphate dihydrate (NaH <sub>2</sub> PO <sub>4</sub> × 2H <sub>2</sub> O):	1.5 mg/ml
Water for injections (H <sub>2</sub> O):	q.s.p. 1.0 ml
The concentration of glycerol in the gel is 17.5 mg/ml.	
Belotero® Revive does not contain any cells, tissues, organs or any other derivatives of animal or human origin.	

#### Indications

Belotero® Revive is indicated to restore skin quality by revitalization of facial skin, characterized by decreased skin hydration, loss of elasticity and firmness, and presence of superficial fine lines.

#### Posology and Administration Method

Belotero® Revive is designed to be injected into mid- to deep dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection to minimize the risk of potential complications.

Inject Belotero® Revive slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles. The risk of an intravascular injection can be reduced by different strategies, including aspiration prior to injection, utilizing lower volumes and serial injections in high-risk areas, treating one side at a time, pinching/tenting the skin to provide more space superficial to the branches of the main arteries, and manual occlusion of the origin of the supratrochlear vessels with the non-dominant finger. Blunt cannulas may reduce, but not eliminate the risk.

General recommended injection techniques are for example: linear or serial threading, fanning, cross-hatching or serial (micro) puncture. The quantity of product to be injected depends on the area to be corrected.

Belotero® Revive must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of Belotero® Revive, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer lock connection during injection and may cause injury to the patient and/or physician. If during the injection, needle disengagement or a gel leak occurs, discard the syringe and the needle and restart the procedure with a new product.

If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle. Never try to straighten a bent needle; instead throw it away and replace it. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct. The graduations on the syringe label are only intended for orientation for the user.

Gently massage the treated area after the injection to distribute the product uniformly.

Before treatment, the patient’s suitability for the treatment and the patient’s need for pain relief (topical anaesthetics, ice packs, distraction techniques, local anaesthetic injections, or nerve

blocks depending on the injection site and size of needle used), should be assessed.

The maximum application volume for Belotero® Revive must not exceed 4 mL per session. Belotero® Revive can be injected under a treatment regimen of 1 to 3 injection sessions. The shortest possible reinjection interval is 4 weeks. This treatment regimen can be renewed at earliest after 6 months, resulting in a maximum annual dose of 24 mL.

#### Contraindications

Belotero® Revive is contra-indicated:

- In case of known hypersensitivity to one of the product’s components, especially to sodium hyaluronate or to glycerol or BDDE,
- In pregnant and breast-feeding women,
- In patients under 18 years old,
- In patients presenting a general infection.
- In patients presenting an active auto-immune disease.
- Do not inject Belotero® Revive into blood vessels. Do not inject Belotero® Revive into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

Do not inject Belotero® Revive into an area previously treated with a permanent dermal filler.

Do not inject Belotero® Revive into the glabellar or nose region.

#### Precautions for Use

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

In the absence of available clinical data on tolerance of the injection of Belotero® Revive in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject Belotero® Revive on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

Belotero® Revive can be used in combination with other Belotero® products during the same session but in different facial areas. Instructions for use of each product should be followed. No clinical data is available on the injection of Belotero® Revive into patient with a Fitzpatrick skin type V/VI.

No clinical data are available on the injection of Belotero® Revive into the hands.

It is recommended not to inject Belotero® Revive in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

Belotero® Revive can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse) only if injected in different facial areas. Practitioners should be experienced, and patients appropriately selected as not only benefits but also adverse events can be cumulative, and causalty of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of Belotero® Revive into an area already treated with other aesthetic products or procedures. Belotero® Revive must not be used in association with other aesthetic techniques such as peeling, dermabrasion or any type of laser treatment before complete healing after the last treatment. In any case, even if the healing occurs earlier, Belotero® Revive must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of Belotero® Revive with the above-mentioned treatments.

Patients using anti-coagulation, anti-platelet, or thrombolytic medications (e.g. warfarin), anti-inflammatory drugs (oral/injectable corticosteroids or non-steroidal anti-inflammatory drugs (NSAIDs; e.g. aspirin, ibuprofen)), or other substances known to increase coagulation time (vitamins or herbal supplements, e.g. Vitamin E, garlic, Ginkgo biloba and St. John’s Wort), from 10 days pre- to 3 days post injection may have increased reactions of hematomas, nodules or bleeding at the injection site.

Injection of Belotero® Revive into patients with a history of previous herpetic eruption may be associated with herpes reactivation and HHV related diseases, (e.g. pityriasis rosea).

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject Belotero® Revive on a case-by-case basis depending of the nature of the disease as well as the associated treatment. Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.

Do not transfer Belotero® Revive into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe. Discard the syringe, the remaining product and the needles in the appropriate container after use.

**Do not re-sterilize and do not reuse due to the associated risks including infection.**

The patient must avoid applying makeup (including skin care products) for at least 12 hours after treatment as well as avoid saunas, peeling, Turkish baths and prolonged exposure to the sun, UV rays, extreme heat and cold for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area and should avoid strenuous physical activity following treatment.

The patient must avoid drinking alcohol for 24 hours before and after treatment. Alcohol may cause the blood vessels to dilate and cause more bruising.

The gel should not be injected into the area where another absorbable implant is located, for example, biodegradable threads for lifting.

#### Warnings

Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore

recommended that Belotero® Revive does not come into contact with such substances.

Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur. Removal or change of the implant is not possible. For accelerated resolution injection of hyaluronidase or longidaze ferments can be used according to the treating physician’s choice. Visualization by MRI and/or ultrasound diagnostics is not possible.

#### Side Effects and Adverse Events

Patients must be informed by the practitioner about possible side effects and adverse events before treatment.

##### Side Effects:

Injection site reactions may occur following injection into the skin but disappear spontaneously within a few days. This includes swelling, nodule or lump/bump, bruising/purpura, hematoma, ecchymosis, induration, erythema/redness, tenderness, pain, discoloration and pruritus/itching, tingling, paresthesia, numbness, hypoesthesia, scabbing, needle mark and discomfort or irritation. These injection site reactions are generally of mild or moderate intensity. A transient bleeding may also occur at the injection site and usually stops spontaneously as soon as the injection is finished.

##### Adverse Events:

In occasional cases, one or more of the following may occur in conjunction with the use of products of the Belotero® portfolio either immediately or as a delayed reaction: acne cystic, milia, skin dryness (rough facial skin, exfoliation), injection site erosion, inflammation, shivering, fatigue, lymphatic system disorder, rash, burning sensation, injection site sore/warmth/fever, pruritus/itching, urticaria, hematoma, telangiectasia, ecchymosis, edema (including lymph edema), headache, cephalgia, tumefaction, tension, swelling (including persistent swelling), hyper- or hypo-pigmentation, angioedema, induration, blister, vesicle, papule, lump/ bump (visible and/or palpable material) or nodule (including inflammatory nodules), mass, granuloma (including inflammatory signs and foreign body reactions), necrosis, ischemia, vascular occlusion, embolization, infarction, Tyndall effect (including translucent chords), hypersensitivity, allergic reactions (including asthma attack, Quincke’s edema, anaphylactic shock or throat tightening) to one of the product’s components (e.g. hyaluronic acid, BDDE), oral and dental disorders, nervous system impairment, impairment of the otorhinolaryngological system (e.g. nasal congestion, oropharyngeal pain, dysgeusia, rhinorrhea, epistaxis, sinusitis, transient hearing loss), mastication pain, parotid gland enlargement, muscle twitching, muscle injury/disorder, nausea, vomiting, circulatory collapse, presyncope, peripheral venous disease, hot flush, anxiety caused by trypanophobia, patient dissatisfaction and disappointment (due to lack of or reduced performance, decreased firmness/response, undesirable aesthetic effect), injection site discharge, device migration, product distribution issue (e.g. product accumulation), injection site indentation, superficial vein prominence, overcorrection or cranial nerve disorder (e.g. cranial nerve paralysis, facial paralysis, trigeminal neuralgia).

Rare cases of the following adverse events have been reported with hyaluronic acid products such as infection (e.g. erysipelas, phlegmon, cellulitis, including open or draining wounds and (dental) abscess, impetigo, pustules) chronic infection (including biofilm formation), scarring, persistent skin discoloration, sensory dysfunction, non-thrombotic lung embolism as well as serous dysfunaoma formation in subjects with hepatic C and interferon treatment, cerebral injuries (e.g. intracranial penetration, subarachnoid hemorrhage), strabismus, ophthalmoplegia, iris adhesions, cataract, conjunctival hemorrhage, eyelid ptosis and laceration.

The risk of granuloma, ischemia, necrosis and vascular occlusion is higher with deep injections and high volumes.

Isolated cases of visual impairment or blindness following unintentional intra-arterial injection have been reported in literature.

Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of colour are more likely to develop post-inflammatory hyperpigmentation and / or hypertrophic scar/keloid formation following injection procedures. Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

There is no known interaction with local or loco-regional anesthetics.

**Circumstances in which the patient should consult a health care professional**

Patients should be instructed to report any side effects which last for more than one week and any adverse event as soon as it occurs to his/her practitioner, especially if patient has changes in his/her vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in his/her face, arms or legs, difficulty walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment. The practitioner may then refer the patient to the appropriate treatment.

#### Assembly of needle to syringe

For optimal use of Belotero® Revive, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

- Firmly hold the glass cylinder of the syringe and the Luer lock adaptor** between the thumb and the index fingers.
- Grasp the protective cap with the other hand and unscrew it.

- Push and twist the needle on the syringe until a resistance is felt.** Do not over-tight. Over-tightening of the needle may lead to the Luer lock moving and dislodging from the syringe.
- Keep holding the Luer lock and remove the sheath from the needle.

#### Transportation, storage and use

Store at temperature between 2 °C to 25 °C. Protect from light and freezing. Avoid mechanical shocks. The manufacturer has no requirements to humidity and atmospheric pressure. It is recommended to perform injection under room temperature (from 15 °C to 25 °C).

The implant is intended for use in the internal environment of the body. The products are explorable during operation when exposed to temperatures from 32 °C to 42 °C. Transportation is carried out in cardboard boxes loaded with consumer packs inside. Shipping by any modes of transport at temperature between 2 °C to 25 °C does not affect product properties.

#### Disposal

In the territory of the Russian Federation used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with SanPIN 2.1.3684-21 „Sanitary and epidemiological requirements to the maintenance of the territory of urban and rural settlements, to waterbodies, to drinking water and water supply of citizens, to atmosphere air, to soils, to living areas, to exploitation of manufacturing and communal areas, to organization and realization of sanitary and anti-epidemic (preventive) events“ as belonging to hazard class B.

#### Requirements for environmental protection

With proper disposal, there is no harm to the environment.